

TCF19

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:20 - 1:100

Calculated MW: 37kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human TCF19

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

SC1; TCF-19;

Catalog #: A6147

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 6941

Isotype: IgG

Swiss Prot: Q9Y242

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein that belongs to the Plant Homeo Domain finger family of transcription factors. The encoded protein is thought to function during the G1/S transition in the cell cycle.

Alternate splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).