

HMGN1

Reactivity:Human Mouse

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW:11kDa

Observed MW:Refer to figures

Immunogen:

A synthetic peptide of human HMGN1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

fq

Synonym:

HMG14;

Catalog #:A6155

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:3150

Isotype:IgG

Swiss Prot:P05114

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene binds nucleosomal DNA and is associated with transcriptionally active chromatin. Along with a similar protein, HMG17, the encoded protein may help maintain an open chromatin configuration around transcribable genes.

To place an order, please [Click HERE](#).