

HMGN2

Reactivity: Human Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 9kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human HMGN2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

f

Synonym:

HMG17;

Catalog #: A6156

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3151

Isotype: IgG

Swiss Prot: P05204

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene binds nucleosomal DNA and is associated with transcriptionally active chromatin. Along with a similar protein, HMGN1, the encoded protein may help maintain an open chromatin configuration around transcribable genes. The protein has also been found to have antimicrobial activity against bacteria, viruses and fungi.

To place an order, please [Click HERE](#).