

IFI6

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 12kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human IFI6

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

6-16; G1P3; FAM14C; IFI616; IFI-6-16;

Catalog #: A6157

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 2537

Isotype: IgG

Swiss Prot: P09912

Purity: Affinity purification

For research use only.

Background:

This gene was first identified as one of the many genes induced by interferon. The encoded protein may play a critical role in the regulation of apoptosis. A minisatellite that consists of 26 repeats of a 12 nucleotide repeating element resembling the mammalian splice donor consensus sequence begins near the end of the second exon. Alternatively spliced transcript variants that encode different isoforms by using the two downstream repeat units as splice donor sites have been described.

To place an order, please [Click HERE](#).