

RNF122

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:17kDa

Observed MW:Refer to Figures

Immunogen:

A synthetic peptide of human RNF122

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

I

Catalog #:A6160

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:79845

Isotype:IgG

Swiss Prot:Q9H9V4

Purity:Affinity purification

For research use only.

Background:

The encoded protein contains a RING finger, a motif present in a variety of functionally distinct proteins and known to be involved in protein-protein and protein-DNA interactions. The encoded protein is localized to the endoplasmic reticulum and golgi apparatus, and may be associated with cell viability.

To place an order, please [Click HERE](#).