

OSM

Reactivity:Human

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW:28kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human OSM

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #:A6163

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:5008

Isotype:IgG

Swiss Prot:P13725

Purity:Affinity purification

For research use only.

Background:

Oncostatin M is a member of a cytokine family that includes leukemia-inhibitory factor, granulocyte colony-stimulating factor, and interleukin 6. This gene encodes a growth regulator which inhibits the proliferation of a number of tumor cell lines. It regulates cytokine production, including IL-6, G-CSF and GM-CSF from endothelial cells.

To place an order, please [Click HERE](#).