

MYLIP

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 50kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human MYLIP

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MIR;MYLIP;

Catalog #: A6166

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 29116

Isotype: IgG

Swiss Prot: Q8WY64

Purity: Affinity purification

For research use only.

Background:

The ERM protein family members ezrin, radixin, and moesin are cytoskeletal effector proteins linking actin to membrane-bound proteins at the cell surface. Myosin regulatory light chain interacting protein (MYLIP) is a novel ERM-like protein that interacts with myosin regulatory light chain and inhibits neurite outgrowth.

To place an order, please [Click HERE](#).