

KLHL25

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 66kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human KLHL25

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ENC2; ENC-2;

To place an order, please [Click HERE](#).

Catalog #: A6167

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 64410

Isotype: IgG

Swiss Prot: Q9H0H3

Purity: Affinity purification

For research use only.