

CSF3

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:1000

Calculated MW: 22kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human CSF3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

GCSF; CSF3OS; C17orf33;

Catalog #: A6178

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1440

Isotype: IgG

Swiss Prot: P09919

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a cytokine that controls the production, differentiation, and function of granulocytes. The active protein is found extracellularly. Alternatively spliced transcript variants have been described for this gene.

To place an order, please [Click HERE](#).