

HDAC10

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 71kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human HDAC10

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HD10;

Catalog #: A6184

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 83933

Isotype: IgG

Swiss Prot: Q969S8

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene belongs to the histone deacetylase family, members of which deacetylate lysine residues on the N-terminal part of the core histones. Histone deacetylation modulates chromatin structure, and plays an important role in transcriptional regulation, cell cycle progression, and developmental events. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).