

CTGF

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:1000 IHC 1:20 - 1:50

Calculated MW: 38kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human CTGF

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

a

Synonym:

CCN2; NOV2; HCS24; IGFBP8

Catalog #: A6190

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1490

Isotype: IgG

Swiss Prot: P29279

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a mitogen that is secreted by vascular endothelial cells. The encoded protein plays a role in chondrocyte proliferation and differentiation, cell adhesion in many cell types, and is related to platelet-derived growth factor. Certain polymorphisms in this gene have been linked with a higher incidence of systemic sclerosis.

To place an order, please [Click HERE](#).