

SMYD5

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF IP

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100 IF 1:50 - 1:100 IP 1:50 - 1:100

Calculated MW: 47kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human SMYD5

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

1 mg/ml

Synonym:

RRG1; RAI15; NN8-4AG; ZMYND23

To place an order, please [Click HERE](#).

Catalog #: A6191

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10322

Isotype: IgG

Swiss Prot: Q6GMV2

Purity: Affinity purification

For research use only.