

TFPI2

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 27kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human TFPI2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PP5; REF1; TFPI-2;

Catalog #: A6214

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 7980

Isotype: IgG

Swiss Prot: P48307

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the Kunitz-type serine proteinase inhibitor family. The protein can inhibit a variety of serine proteases including factor VIIa/tissue factor, factor Xa, plasmin, trypsin, chymotrypsin and plasma kallikrein. This gene has been identified as a tumor suppressor gene in several types of cancer. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).