

FLOT1

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:10 - 1:100

Calculated MW: 47kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human FLOT1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A6220

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10211

Isotype: IgG

Swiss Prot: O75955

Purity: Affinity purification

For research use only.

Background:

Caveolae are small domains on the inner cell membrane involved in vesicular trafficking and signal transduction. FLOT1 encodes a caveolae-associated, integral membrane protein. The function of flotillin 1 has not been determined.

To place an order, please [Click HERE](#).