

CBX8

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:10 - 1:100

Calculated MW: 43kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human CBX8

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

fl

Synonym:

PC3; RC1;

To place an order, please [Click HERE](#).

Catalog #: A6222

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 57332

Isotype: IgG

Swiss Prot: Q9HC52

Purity: Affinity purification

For research use only.