

CD5L

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 38kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human CD5L

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

AIM; API6; PRO229; Spalpha; SP-ALPHA;

To place an order, please [Click HERE](#).

Catalog #: A6223

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 922

Isotype: IgG

Swiss Prot: O43866

Purity: Affinity purification

For research use only.