

GIP

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 17kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human GIP

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A6230

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 2695

Isotype: IgG

Swiss Prot: P09681

Purity: Affinity purification

For research use only.

Background:

This gene encodes an incretin hormone and belongs to the glucagon superfamily. The encoded protein is important in maintaining glucose homeostasis as it is a potent stimulator of insulin secretion from pancreatic beta-cells following food ingestion and nutrient absorption. This gene stimulates insulin secretion via its G protein-coupled receptor activation of adenylyl cyclase and other signal transduction pathways. It is a relatively poor inhibitor of gastric acid secretion.

To place an order, please [Click HERE](#).