

NR2E3

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 44kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human NR2E3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PNR; RNR; rd7; ESCS; RP37;

Catalog #: A6234

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10002

Isotype: IgG

Swiss Prot: Q9Y5X4

Purity: Affinity purification

For research use only.

Background:

This protein is part of a large family of nuclear receptor transcription factors involved in signaling pathways. Nuclear receptors have been shown to regulate pathways involved in embryonic development, as well as in maintenance of proper cell function in adults. Members of this family are characterized by discrete domains that function in DNA and ligand binding. This gene encodes a retinal nuclear receptor that is a ligand-dependent transcription factor. Defects in this gene are a cause of enhanced S cone syndrome. Alternatively spliced transcript variants encoding different isoforms have been identified.

To place an order, please [Click HERE](#).