

AKT1S1

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 27kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human AKT1S1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

Lobe; PRAS40;

Catalog #: A6238

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 84335

Isotype: IgG

Swiss Prot: Q96B36

Purity: Affinity purification

For research use only.

Background:

AKT1S1 is a proline-rich substrate of AKT (MIM 164730) that binds 14-3-3 protein (see YWHAH, MIM 113508) when phosphorylated (Kovacina et al., 2003)

To place an order, please [Click HERE](#).