

PHLDA2

Reactivity: Human

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:20 - 1:100

Calculated MW: 17kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human PHLDA2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

IPL; BRW1C; BWR1C; HLDA2; TSSC3;

Catalog #: A6244

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 7262

Isotype: IgG

Swiss Prot: Q53GA4

Purity: Affinity purification

For research use only.

Background:

This gene is located in a cluster of imprinted genes on chromosome 11p15.5, which is considered to be an important tumor suppressor gene region. Alterations in this region may be associated with the Beckwith-Wiedemann syndrome, Wilms tumor, rhabdomyosarcoma, adrenocortical carcinoma, and lung, ovarian, and breast cancer. This gene has been shown to be imprinted, with preferential expression from the maternal allele in placenta and liver.

To place an order, please [Click HERE](#).