

PNPLA2

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 55kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human PNPLA2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ATGL; TTS2; PEDF-R; FP17548; TTS-2.2; iPLA2zeta; 1110001C14Rik;

Catalog #: A6245

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 57104

Isotype: IgG

Swiss Prot: Q96AD5

Purity: Affinity purification

For research use only.

Background:

This gene encodes an enzyme which catalyzes the first step in the hydrolysis of triglycerides in adipose tissue. Mutations in this gene are associated with neutral lipid storage disease with myopathy.

To place an order, please [Click HERE](#).