

AOC1

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 85kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human AOC1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ABP; DAO; KAO; ABP1; DAO1;

Catalog #: A6249

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 26

Isotype: IgG

Swiss Prot: P19801

Purity: Affinity purification

For research use only.

Background:

This gene encodes a metal-binding membrane glycoprotein that oxidatively deaminates putrescine, histamine, and related compounds. The encoded protein is inhibited by amiloride, a diuretic that acts by closing epithelial sodium ion channels. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.

To place an order, please [Click HERE](#).