

CLEC12A

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 30kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human CLEC12A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CLL1; MICL; CLL-1; DCAL-2;

Catalog #: A6250

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 160364

Isotype: IgG

Swiss Prot: Q5QGZ9

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the C-type lectin/C-type lectin-like domain (CTL/CTLD) superfamily. Members of this family share a common protein fold and have diverse functions, such as cell adhesion, cell-cell signaling, glycoprotein turnover, and roles in inflammation and immune response. The protein encoded by this gene is a negative regulator of granulocyte and monocyte function. Several alternatively spliced transcript variants of this gene have been described, but the full-length nature of some of these variants has not been determined. This gene is closely linked to other CTL/CTLD superfamily members in the natural killer gene complex region on chromosome 12p13.

To place an order, please [Click HERE](#).