

SPRY4

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 33kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human SPRY4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

q

Synonym:

HH17;

Catalog #: A6263

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 81848

Isotype: IgG

Swiss Prot: Q9C004

Purity: Affinity purification

For research use only.

Background:

SPRY4 is an inhibitor of the receptor-transduced mitogen-activated protein kinase (MAPK) signaling pathway. It is positioned upstream of RAS (see HRAS; MIM 190020) activation and impairs the formation of active GTP-RAS.

To place an order, please [Click HERE](#).