

TNFRSF10A

Reactivity: Human Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:200 - 1:500 IHC 1:50 - 1:100

Calculated MW: 50kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human TNFRSF10A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

bp

Synonym:

DR4; APO2; CD261; TRAILR1; TRAILR-1

Catalog #: A6267

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 8797

Isotype: IgG

Swiss Prot: O00220

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of the TNF-receptor superfamily. This receptor is activated by tumor necrosis factor-related apoptosis inducing ligand (TNFSF10/TRAIL), and thus transduces cell death signal and induces cell apoptosis. Studies with FADD-deficient mice suggested that FADD, a death domain containing adaptor protein, is required for the apoptosis mediated by this protein.

To place an order, please [Click HERE](#).