

LRP12

Reactivity: Human

Tested applications: WB IHC IF FC

Calculated MW: 859

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human LRP12

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #: A6275

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 29967

Isotype: IgG

Swiss Prot: Q9Y561

Purity: Affinity purification

For research use only.