

AGFG1

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 58kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human AGFG1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HRB; RAB; RIP;

Catalog #: A6294

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3267

Isotype: IgG

Swiss Prot: P52594

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is related to nucleoporins, a class of proteins that mediate nucleocytoplasmic transport. The encoded protein binds the activation domain of the human immunodeficiency virus Rev protein when Rev is assembled onto its RNA target, and is required for the nuclear export of Rev-directed RNAs. Multiple transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).