

OXA1L

Reactivity: Human Mouse

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:10 - 1:100

Calculated MW: 48kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human OXA1L

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

OXA1;

To place an order, please [Click HERE](#).

Catalog #: A6300

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 5018

Isotype: IgG

Swiss Prot: Q15070

Purity: Affinity purification

For research use only.