

SIVA1

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 19kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human SIVA1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

SIVA; CD27BP; Siva-1; Siva-2;

Catalog #: A6326

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10572

Isotype: IgG

Swiss Prot: O15304

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein with an important role in the apoptotic (programmed cell death) pathway induced by the CD27 antigen, a member of the tumor necrosis factor receptor (TNFR) superfamily. The CD27 antigen cytoplasmic tail binds to the N-terminus of this protein. Two alternatively spliced transcript variants encoding distinct proteins have been described.

To place an order, please [Click HERE](#).