

## INTS4

**Reactivity:** Human Mouse Rat

**Tested applications:** WB IHC

**Recommended Dilution:** WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:** 108kDa

**Observed MW:** Refer to Figures

**Immunogen:**

Recombinant protein of human INTS4

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

INT4; MST093;

**Catalog #:** A6339

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 92105

**Isotype:** IgG

**Swiss Prot:** Q96HW7

**Purity:** Affinity purification

For research use only.

**Background:**

INTS4 is a subunit of the Integrator complex, which associates with the C-terminal domain of RNA polymerase II large subunit (POLR2A; MIM 180660) and mediates 3-prime end processing of small nuclear RNAs U1 (RNU1; MIM 180680) and U2 (RNU2; MIM 180690) (Baillat et al., 2005)

*To place an order, please [Click HERE](#).*