

SYP

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF IP

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100 IP 1:20 - 1:50

Calculated MW: 34kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human SYP

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MRX96; MRXSYP;

Catalog #: A6344

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 6855

Isotype: IgG

Swiss Prot: P08247

Purity: Affinity purification

For research use only.

Background:

This gene encodes an integral membrane protein of small synaptic vesicles in brain and endocrine cells. The protein also binds cholesterol and is thought to direct targeting of vesicle-associated membrane protein 2 (synaptobrevin) to intracellular compartments. Mutations in this gene are associated with X-linked mental retardation (XLMR).

To place an order, please [Click HERE](#).