

C1R

Reactivity:Human Mouse

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW:80kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human C1R

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #:A6360

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:715

Isotype:IgG

Swiss Prot:P00736

Purity:Affinity purification

For research use only.

To place an order, please [Click HERE](#).