

CLIC1

Reactivity:Human

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:10 - 1:100

Calculated MW:26kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human CLIC1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

G6; NCC27;

Catalog #:A6363

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:1192

Isotype:IgG

Swiss Prot:O00299

Purity:Affinity purification

For research use only.

Background:

Chloride channels are a diverse group of proteins that regulate fundamental cellular processes including stabilization of cell membrane potential, transepithelial transport, maintenance of intracellular pH, and regulation of cell volume. Chloride intracellular channel 1 is a member of the p64 family; the protein localizes principally to the cell nucleus and exhibits both nuclear and plasma membrane chloride ion channel activity.

To place an order, please [Click HERE](#).