

## CLPS

**Reactivity:** Human Mouse

**Tested applications:** WB IHC IF

**Recommended Dilution:** WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:10 - 1:100

**Calculated MW:** 11kDa

**Observed MW:** Refer to Figures

**Immunogen:**

Recombinant protein of human CLPS

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Catalog #:** A6364

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 1208

**Isotype:** IgG

**Swiss Prot:** P04118

**Purity:** Affinity purification

For research use only.

**Background:**

The protein encoded by this gene is a cofactor needed by pancreatic lipase for efficient dietary lipid hydrolysis. It binds to the C-terminal, non-catalytic domain of lipase, thereby stabilizing an active conformation and considerably increasing the overall hydrophobic binding site. The gene product allows lipase to anchor noncovalently to the surface of lipid micelles, counteracting the destabilizing influence of intestinal bile salts. This cofactor is only expressed in pancreatic acinar cells, suggesting regulation of expression by tissue-specific elements. Three transcript variants encoding different isoforms have been found for this gene.

*To place an order, please [Click HERE](#).*