

SLC26A2

Reactivity: Human

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 81kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human SLC26A2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

DTD; EDM4; DTDST; MST153; D5S1708; MSTP157;

Catalog #: A6369

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1836

Isotype: IgG

Swiss Prot: P50443

Purity: Affinity purification

For research use only.

Background:

The diastrophic dysplasia sulfate transporter is a transmembrane glycoprotein implicated in the pathogenesis of several human chondrodysplasias. It apparently is critical in cartilage for sulfation of proteoglycans and matrix organization.

To place an order, please [Click HERE](#).