

ELF3

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:10 - 1:100

Calculated MW: 41kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human ELF3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ERT; ESX; EPR-1; ESE-1;

To place an order, please [Click HERE](#).

Catalog #: A6371

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1999

Isotype: IgG

Swiss Prot: P78545

Purity: Affinity purification

For research use only.