

GRK6

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:10 - 1:100

Calculated MW: 65kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human GRK6

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

GPRK6;

Catalog #: A6379

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 2870

Isotype: IgG

Swiss Prot: P43250

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the guanine nucleotide-binding protein (G protein)-coupled receptor kinase subfamily of the Ser/Thr protein kinase family. The protein phosphorylates the activated forms of G protein-coupled receptors thus initiating their deactivation. Several transcript variants encoding different isoforms have been described for this gene.

To place an order, please [Click HERE](#).