

LIPA

Reactivity:Human

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:45kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human LIPA

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

LAL; CESD;

Catalog #:A6385

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:3988

Isotype:IgG

Swiss Prot:P38571

Purity:Affinity purification

For research use only.

Background:

This gene encodes lipase A, the lysosomal acid lipase (also known as cholesterol ester hydrolase). This enzyme functions in the lysosome to catalyze the hydrolysis of cholesteryl esters and triglycerides. Mutations in this gene can result in Wolman disease and cholesteryl ester storage disease. Alternatively spliced transcript variants have been found for this gene.

To place an order, please [Click HERE](#).