

MLN

Reactivity: Human Mouse

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100 IF 1:50 - 1:100

Calculated MW: 12kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human MLN

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A6388

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4295

Isotype: IgG

Swiss Prot: P12872

Purity: Affinity purification

For research use only.

Background:

This gene encodes a small peptide hormone that is secreted by cells of the small intestine to regulate gastrointestinal contractions and motility. Proteolytic processing of the secreted protein produces the mature peptide and a byproduct referred to as motilin-associated peptide (MAP). Three transcript variants encoding different preproprotein isoforms but the same mature peptide have been found for this gene.

To place an order, please [Click HERE](#).