

MSRA

Reactivity: Human

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:10 - 1:100

Calculated MW: 26kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human MSRA

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PMSR;

Catalog #: A6389

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4482

Isotype: IgG

Swiss Prot: Q9UJ68

Purity: Affinity purification

For research use only.

Background:

This gene encodes a ubiquitous and highly conserved protein that carries out the enzymatic reduction of methionine sulfoxide to methionine. Human and animal studies have shown the highest levels of expression in kidney and nervous tissue. The protein functions in the repair of oxidatively damaged proteins to restore biological activity. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).