

PNLIP

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 51kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human PNLIP

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PL; PTL; PNLIPD;

Catalog #: A6396

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 5406

Isotype: IgG

Swiss Prot: P16233

Purity: Affinity purification

For research use only.

Background:

This gene is a member of the lipase gene family. It encodes a carboxyl esterase that hydrolyzes insoluble, emulsified triglycerides, and is essential for the efficient digestion of dietary fats. This gene is expressed specifically in the pancreas.

To place an order, please [Click HERE](#).