

TRIM27

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:10 - 1:100

Calculated MW: 58kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human TRIM27

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

RFP; RNF76;

Catalog #: A6405

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 5987

Isotype: IgG

Swiss Prot: P14373

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the tripartite motif (TRIM) family. The TRIM motif includes three zinc-binding domains, a RING, a B-box type 1 and a B-box type 2, and a coiled-coil region. This protein localizes to the nuclear matrix. It interacts with the enhancer of polycomb protein and represses gene transcription. It is also thought to be involved in the differentiation of male germ cells. Fusion of the N-terminus of this protein with the truncated C-terminus of the RET gene product has been shown to result in production of the ret transforming protein.

To place an order, please [Click HERE](#).