

XCL1

Reactivity: Human

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 12kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human XCL1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

LTN; ATAC; LPTN; SCM1; SCM-1; SCM1A; SCYC1; SCM-1a;

Catalog #: A6408

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 6375

Isotype: IgG

Swiss Prot: P47992

Purity: Affinity purification

For research use only.

Background:

This antimicrobial gene encodes a member of the chemokine superfamily. Chemokines function in inflammatory and immunological responses, inducing leukocyte migration and activation. The encoded protein is a member of the C-chemokine subfamily, retaining only two of four cysteines conserved in other chemokines, and is thought to be specifically chemotactic for T cells. This gene and a closely related family member are located on the long arm of chromosome 1.

To place an order, please [Click HERE](#).