

SNRPA

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 31kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human SNRPA

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

U1A; Mud1; U1-A;

Catalog #: A6410

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 6626

Isotype: IgG

Swiss Prot: P09012

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene associates with stem loop II of the U1 small nuclear ribonucleoprotein, which binds the 5' splice site of precursor mRNAs and is required for splicing.

The encoded protein autoregulates itself by polyadenylation inhibition of its own pre-mRNA via dimerization and has been implicated in the coupling of splicing and polyadenylation.

To place an order, please [Click HERE](#).