

SPINT1

Reactivity: Human Mouse

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100 IF 1:50 - 1:100

Calculated MW: 58kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human SPINT1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HAI; HAI1; MANSC2;

Catalog #: A6411

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 6692

Isotype: IgG

Swiss Prot: O43278

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of the Kunitz family of serine protease inhibitors.

The protein is a potent inhibitor specific for HGF activator and is thought to be involved in the regulation of the proteolytic activation of HGF in injured tissues. Alternative splicing results in multiple variants encoding different isoforms.

To place an order, please [Click HERE](#).