

SPN

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 40kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human SPN

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

LSN; CD43; GALGP; GPL115;

Catalog #: A6412

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 6693

Isotype: IgG

Swiss Prot: P16150

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a major sialoglycoprotein found on the surface of thymocytes, T lymphocytes, monocytes, granulocytes, and some B lymphocytes. It may be part of a physiologic ligand-receptor complex involved in T-cell activation. During T-cell activation, this protein is actively removed from the T-cell-APC (antigen-presenting cell) contact site, suggesting a negative regulatory role in adaptive immune response.

To place an order, please [Click HERE](#).