

KLF10

Reactivity: Human

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 52kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human KLF10

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

EGRA; TIEG; TIEG1; EGR-alpha;

Catalog #: A6415

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 7071

Isotype: IgG

Swiss Prot: Q13118

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of a family of proteins that feature C2H2-type zinc finger domains. The encoded protein is a transcriptional repressor that acts as an effector of transforming growth factor beta signaling. Activity of this protein may inhibit the growth of cancers, particularly pancreatic cancer. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).