

TIMP4

Reactivity: Human

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:10 - 1:100

Calculated MW: 25kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human TIMP4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A6416

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 7079

Isotype: IgG

Swiss Prot: Q99727

Purity: Affinity purification

For research use only.

Background:

This gene belongs to the TIMP gene family. The proteins encoded by this gene family are inhibitors of the matrix metalloproteinases, a group of peptidases involved in degradation of the extracellular matrix. The secreted, netrin domain-containing protein encoded by this gene is involved in regulation of platelet aggregation and recruitment and may play role in hormonal regulation and endometrial tissue remodeling.

To place an order, please [Click HERE](#).