

## CLEC3B

---

**Reactivity:** Human

**Tested applications:** WB IHC IF

**Recommended Dilution:** WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:10 - 1:100

**Calculated MW:** 22kDa

**Observed MW:** Refer to Figures

**Immunogen:**

Recombinant protein of human CLEC3B

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

TN; TNA;

*To place an order, please [Click HERE](#).*

**Catalog #:** A6418

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 7123

**Isotype:** IgG

**Swiss Prot:** P05452

**Purity:** Affinity purification

For research use only.