

SF1

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 68kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human SF1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

BBP; MBBP; ZFM1; ZNF162; D11S636; ZCCHC25;

Catalog #: A6424

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 7536

Isotype: IgG

Swiss Prot: Q15637

Purity: Affinity purification

For research use only.

Background:

This gene encodes a nuclear pre-mRNA splicing factor. The encoded protein specifically recognizes the intron branch point sequence and is required for the early stages of spliceosome assembly. Alternate splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).